



**EasyThru<sup>®</sup>**

FOR MEDICAL TECHNOLOGY  
FOR LIFE CARE

**2021 version**

**YILI MEDICAL**



[www.yili-medical.com](http://www.yili-medical.com)



0197

A modern, multi-story building with a grey metallic facade and large windows. The building is identified by a large logo on its upper right side. A thick green curved line sweeps across the left side of the image, separating the header area from the main content.

**YILI  
Medical**

# ABOUT US

**YILI MEDICAL** is a professional manufacturer of medical disposables, established in 2010

**YILI MEDICAL** is a professional manufacturer and supplier of medical products over 10 years. With the rapid development of sales and the demand of the expanding market, YiLi medical company has grown from a single product to diversified products. Our main products are including breathing anesthesia product, urinary, and injection infusion etc. Three family owned factories are running different product lines ensure we are able to offer various products and offer you one-stop service.

# HISTORY OF YILI MEDICAL

- 2010 Establish Medical Research and Development Institution
- 2011 Establish YiLi Medical Sales Department in Hong Kong
- 2012 Establish Nanchang YiLi Medical Instrument Co., Ltd in Nanchang City
- 2013 Establish first factory in Jinxian town of Nanchang
- 2014 got ISO and CE certificated, the same year dispatched sales team to Mexico and Turkey for marketing survey
- 2016 Invest other two family own factories to expand products range and also registered EasyThru Brand
- 2017 Sign distributors in Turkey, Korea, Germany
- 2018 FDA certificate for Foley catheter and Oxygen Mask, and become a sub sidiary of YiKang Medical Group.
- 2019 Win tenders in Qatar, UK
- 2020 Registered disposable face mask in Germany

WE HAVE WORKED IN MEDICAL FIELD MORE THAN 12 YEARS



2014



2016



2015



2017

Here we are



## PVC Products Production Line

### FACTORY SITE INFORMATION

- Factory covers 80000 square meters
- 100000 level cleaning room 80000 square meters
- Workers 1200

### QUALITY CONTROL

- Production condition is 100000 cleaning room
- Factories are ISO 13485 and CE certificated
- Each process is inspected and trackable, from raw material coming, semi products, packaging, testing, packing, sterilization, loading container, after sales service. We focus on every details.



1: Tube Extrusion



2: Printing



4: Assembling



3: Balloon Testing

## Latex Foley Catheter Production Line



1: Imported raw material



2: Auto - dipping



3: Eyelets punching



4: Clean



8: Packing



7: ETesting Balloon




6: Assembling



5: Measure size

# CERTIFIED PRODUCT

  
TÜVRheinland

**EC Certificate**  
Directive 93/42/EEC Annex II, excluding Section 4  
Full Quality Assurance System  
Medical Devices

Registration No.: HD 60146269 0001  
Report No.: 50307115 002


**Manufacturer:** Nanchang Yili Medical Instrument Co., Ltd.  
Room 908, 101# Building, JLH605-B01  
Greenland International Exp Center  
Jiulong Road 1177, Honggutan District  
330000 Nanchang, Jiangxi  
P.R. China

**Products:** Disposable Urine Catheter for Single Use

**Expiry Date:** 2024-05-20

The Notified Body hereby declares that the requirements of Annex II, excluding section 4 of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex II, section 5 of the aforementioned directive. For placing on the market of class III devices covered by this certificate an EC design-examination certificate according to Annex II, section 4 is required.

**Effective Date:** 2020-04-13  
**Date:** 2020-04-13

  
Notified Body:  
Wenxiang Zhang

**TÜV Rheinland LGA Products GmbH · Tillystraße 2 · 90431 Nürnberg**  
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

  
TÜVRheinland

**Certificate**  
The Certification Body of  
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization  
**Nanchang Yili Medical Instrument Co., Ltd**  
Room 908, 101# Building, JLH605-B01  
Greenland International Exp Center  
Jiulong Road 1177, Honggutan District  
330000 Nanchang, Jiangxi  
P.R. China

has established and applies a quality management system for medical devices for the following scope  
(See attachment for Scope included)

Proof has been furnished that the requirements specified in  
**EN ISO 13485:2016**  
are fulfilled. The quality management system is subject to yearly surveillance

**Effective Date:** 2020-04-13  
**Certificate Registration No.:** SX 60146270 0001  
**An audit was performed Report No.:** 50307115 002  
**This Certificate is valid until:** 2023-02-12

Certification Body

  
DAKKS  
Dachdecker  
Ausschuss  
D 24 14369 21 02

  
Wenxiang Zhang  
Notified Body

**Date:** 2020-04-13

**TÜV Rheinland LGA Products GmbH · Tillystraße 2 · 90431 Nürnberg**  
Tel: +49 91 866 1371 Fax: +49 91 866 1363 E-Mail: info@tuev.com Web: www.tuev.com



**EC Certificate**  
Directive 93/42/EEC Annex V  
Production Quality Assurance  
Medical Devices

Registration No.: DD 60146267 0001

Report No.: 50307115 002

**Manufacturer:** Nanchang Yili Medical Instrument  
Co., Ltd  
Room 908, 101# Building, JLH605-B01  
Greenland International Exp Center  
Jiulong Road 1177, Honggutan District  
330000 Nanchang, Jiangxi  
P.R. China

**Products:** Medical Devices  
(See Attachment for Products included)

**Expiry Date:** 2024-05-26

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type examination certificate according to Annex III is required.

**Effective Date:** 2020-04-13

**Date:** 2020-04-13

Notified Body



Wenxiang Zhang

**TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg**  
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC  
concerning medical devices with the identification number 0197.



**TÜV Rheinland**  
**LGA Products GmbH**  
Tillystraße 2, 90431 Nürnberg

Doc 17.1.1 Rev. 0

**Attachment to**  
**Certificate**  
**Registration No.:** SX 60146270 0001  
**Report No.:** 50307115 002

**Organization:** Nanchang Yili Medical Instrument  
Co., Ltd  
Room 908, 101# Building, JLH605-B01  
Greenland International Exp Center  
Jiulong Road 1177, Honggutan District  
330000 Nanchang, Jiangxi  
P.R. China

**Scope:** Scope:

Design and Development, Manufacture and Distribution of  
Syringe Needles, Catheters for Single Use,  
Manufacture and Distribution of Insulin Pen Needles,  
Closed Suction Catheters, Arterial Tracheal Tubes for  
Single Use, Disposable Silicone Intubation Masks Disposable  
Circumferential Sheaths, Suction Catheters, Urine Bags and  
VAC External Tubes

**Certification Body**



**Date:** 2020-04-13



Wenxiang Zhang

# PRODUCT LIST

**EasyThru®**

## ANESTHESIA

Endotracheal Tube Cuffed  
Endotracheal Tube Uncuffed  
Oral Preformed Endotracheal Tube Plain  
Oral Preformed Endotracheal Tube CUFF  
Nasal Preformed Endotracheal Tube Plain  
Nasal Preformed Endotracheal Tube Cuffed  
Reinforced Endotracheal Tube Cuffed  
Reinforced endotracheal tube uncuffed  
Endotracheal Tube with Suction Lumen  
Silicone Laryngeal Mask  
Reusable / Reinforced Laryngeal Mask  
PVC Laryngeal Mask  
Enhanced Suction Tube  
Disposable Closed Suction Catheter  
PVC suction catheter  
Suction Connecting Tube + Yankauer (Suction Set)  
PVC Ryle Tube / Nasogastric Tube  
Tracheostomy Tube with Cuff  
Tracheostomy Tube without Cuff  
Oxygen Mask  
Tracheostomy Mask

## OTHERS

N95 Dusty Mask  
Disposable Medical Face Mask (TYPE IIR, Earloop) Surgical Face Mask  
Disposable Medical Face Mask (TYPE IIR, Tie) Surgical Face Mask  
Disposable Surgical Gloves

## URETHRAL PRODUCTS

2 Way Latex Foley Catheter with Silicone Coated  
3 Way Latex Foley Catheter with Silicone Coated  
Haemostatic Catheter (Couvelaire tip Dufour Tip , Mercier tip)  
Double Balloon 3 Way Latex Urinary Catheter  
2 Way Tiemann Tip Latex Foley Catheter  
2 Way Silicone Foley Catheter  
3 Way Silicone Foley Catheter  
PVC Intermittent Urethral catheter(PVC Nelaton catheter)  
Fire-Polished Intermittent Urethral Catheter  
Plastibell Circumcision  
Circumcision Device  
Urine Drainage Bag  
Male External Catheter

## INJECTION

Insulin Pen Needle  
Insulin Syringe  
Safety Blood Lancet  
Needle Holder

## ENTERAL FEEDING PRODUCTS

Silicone Stomach Tube  
TPU Feeding Tube  
PVC Feeding Tube



ENDOTRACHEAL TUBE CUFFED



- Features:
- 1. Medical grade PVC, Non-toxic, Latex free
  - 2. X- ray through the tube
  - 3. Smooth Murphy Eyelet
  - 4. High volume low pressure balloon

01

Product NO	Cuff Size	Box/CTN	Product NO	Cuff Size	Box/CTN
YL-ETTC-20	2.0	10pc/box, 100pcs/CTN	YL-ETTC-65	6.5	10pc/box, 100pcs/CTN
YL-ETTC-25	2.5	10pc/box, 100pcs/CTN	YL-ETTC-70	7.0	10pc/box, 100pcs/CTN
YL-ETTC-30	3.0	10pc/box, 100pcs/CTN	YL-ETTC-75	7.5	10pc/box, 100pcs/CTN
YL-ETTC-35	3.5	10pc/box, 100pcs/CTN	YL-ETTC-80	8.0	10pc/box, 100pcs/CTN
YL-ETTC-40	4.0	10pc/box, 100pcs/CTN	YL-ETTC-85	8.5	10pc/box, 100pcs/CTN
YL-ETTC-45	4.5	10pc/box, 100pcs/CTN	YL-ETTC-90	9.0	10pc/box, 100pcs/CTN
YL-ETTC-50	5.0	10pc/box, 100pcs/CTN	YL-ETTC-95	9.5	10pc/box, 100pcs/CTN
YL-ETTC-55	5.5	10pc/box, 100pcs/CTN	YL-ETTC-100	10.0	10pc/box, 100pcs/CTN
YL-ETTC-60	6.0	10pc/box, 100pcs/CTN	YL-ETTC-105	10.5	10pc/box, 100pcs/CTN

ENDOTRACHEAL TUBE UNCUFFED



- Features:
- 1. Medical grade PVC, Non-toxic, Latex free
  - 2. X- ray through the tube
  - 3. Smooth Murphy Eyelet
  - 4. High volume low pressure balloon
  - 5. Without Balloon

02

Product NO	Cuff Size	Box/CTN	Product NO	Cuff Size	Box/CTN
YL-ETTU-20	2.0	10pc/box, 100pcs/CTN	YL-ETTU-65	6.5	10pc/box, 100pcs/CTN
YL-ETTU-25	2.5	10pc/box, 100pcs/CTN	YL-ETTU-70	7.0	10pc/box, 100pcs/CTN
YL-ETTU-30	3.0	10pc/box, 100pcs/CTN	YL-ETTU-75	7.5	10pc/box, 100pcs/CTN
YL-ETTU-35	3.5	10pc/box, 100pcs/CTN	YL-ETTU-80	8.0	10pc/box, 100pcs/CTN
YL-ETTU-40	4.0	10pc/box, 100pcs/CTN	YL-ETTU-85	8.5	10pc/box, 100pcs/CTN
YL-ETTU-45	4.5	10pc/box, 100pcs/CTN	YL-ETTU-90	9.0	10pc/box, 100pcs/CTN
YL-ETTU-50	5.0	10pc/box, 100pcs/CTN	YL-ETTU-95	9.5	10pc/box, 100pcs/CTN
YL-ETTU-55	5.5	10pc/box, 100pcs/CTN	YL-ETTU-100	10.0	10pc/box, 100pcs/CTN
YL-ETTU-60	6.0	10pc/box, 100pcs/CTN	YL-ETTU-105	10.5	10pc/box, 100pcs/CTN

ORAL PREFORMED ENDOTRACHEAL TUBE PLAIN



- Features :
- 1. Special designed tube shape for oral intubation, anti-kink PVC material
  - 2. Medical grade PVC, Non-toxic, Latex free
  - 3. X- ray through the tube
  - 4. Smooth Murphy Eyelet
  - 5. Without Balloon

Product NO	Cuff Size	Box/CTN	Product NO	Cuff Size	Box/CTN
YL-ETTUAO-30	3.0	10pc/box, 100pcs/CTN	YL-ETTUAO-65	6.5	10pc/box, 100pcs/CTN
YL-ETTUAO-35	3.5	10pc/box, 100pcs/CTN	YL-ETTUAO-70	7.0	10pc/box, 100pcs/CTN
YL-ETTUAO-40	4.0	10pc/box, 100pcs/CTN	YL-ETTUAO-75	7.5	10pc/box, 100pcs/CTN
YL-ETTUAO-45	4.5	10pc/box, 100pcs/CTN	YL-ETTUAO-80	8.0	10pc/box, 100pcs/CTN
YL-ETTUAO-50	5.0	10pc/box, 100pcs/CTN	YL-ETTUAO-85	8.5	10pc/box, 100pcs/CTN
YL-ETTUAO-55	5.5	10pc/box, 100pcs/CTN	YL-ETTUAO-90	9.0	10pc/box, 100pcs/CTN
YL-ETTUAO-60	6.0	10pc/box, 100pcs/CTN	YL-ETTUAO-95	9.5	10pc/box, 100pcs/CTN

ORAL PREFORMED ENDOTRACHEAL TUBE CUFF



- Features :
- 1. Special designed tube shape for oral intubation, anti-kink PVC material
  - 2. Medical grade PVC, Non-toxic, Latex free
  - 3. X- ray through the tube
  - 4. Smooth Murphy Eyelet
  - 5. With Balloon

04

Product NO	Cuff Size	Box/CTN	Product NO	Cuff Size	Box/CTN
YL-ETTCO-30	3.0	10pc/box, 100pcs/CTN	YL-ETTCO-65	6.5	10pc/box, 100pcs/CTN
YL-ETTCO-35	3.5	10pc/box, 100pcs/CTN	YL-ETTCO-70	7.0	10pc/box, 100pcs/CTN
YL-ETTCO-40	4.0	10pc/box, 100pcs/CTN	YL-ETTCO-75	7.5	10pc/box, 100pcs/CTN
YL-ETTCO-45	4.5	10pc/box, 100pcs/CTN	YL-ETTCO-80	8.0	10pc/box, 100pcs/CTN
YL-ETTCO-50	5.0	10pc/box, 100pcs/CTN	YL-ETTCO-85	8.5	10pc/box, 100pcs/CTN
YL-ETTCO-55	5.5	10pc/box, 100pcs/CTN	YL-ETTCO-90	9.0	10pc/box, 100pcs/CTN
YL-ETTCO-60	6.0	10pc/box, 100pcs/CTN	YL-ETTCO-95	9.5	10pc/box, 100pcs/CTN

NASAL PREFORMED ENDOTRACHEAL TUBE PLAIN



- Features :
- 1. Special designed tube shape for Nasal intubation, anti-kink PVC material
  - 2. X- ray through the tube
  - 3. Smooth Murphy Eyelet
  - 4. Without Balloon

Product NO	Cuff Size	Box/CTN	Product NO	Cuff Size	Box/CTN
YL-ETTUN -20	2.0	10pc/box, 100pcs/CTN	YL-ETTUN-65	6.5	10pc/box, 100pcs/CTN
YL-ETTUN-25	2.5	10pc/box, 100pcs/CTN	YL-ETTUN-70	7.0	10pc/box, 100pcs/CTN
YL-ETTUN-30	3.0	10pc/box, 100pcs/CTN	YL-ETTUN-75	7.5	10pc/box, 100pcs/CTN
YL-ETTUN-35	3.5	10pc/box, 100pcs/CTN	YL-ETTUN-80	8.0	10pc/box, 100pcs/CTN
YL-ETTUN-40	4.0	10pc/box, 100pcs/CTN	YL-ETTUN-85	8.5	10pc/box, 100pcs/CTN
YL-ETTUN-45	4.5	10pc/box, 100pcs/CTN	YL-ETTUN-90	9.0	10pc/box, 100pcs/CTN
YL-ETTUN-50	5.0	10pc/box, 100pcs/CTN	YL-ETTUN-95	9.5	10pc/box, 100pcs/CTN
YL-ETTUN-55	5.5	10pc/box, 100pcs/CTN	YL-ETTUN-100	10.0	10pc/box, 100pcs/CTN
YL-ETTUN-60	6.0	10pc/box, 100pcs/CTN			



NASAL PREFORMED ENDOTRACHEAL TUBE (CUFFED)



- Features :
- 1. Special designed tube shape for Nasal intubation, anti-kink PVC material
  - 2. X- ray through the tube
  - 3. Smooth Murphy Eyelet
  - 4. With balloon

06

Product NO	Cuff Size	Box/CTN	Product NO	Cuff Size	Box/CTN
YL-ETTCN-30	3.0	10pc/box, 100pcs/CTN	YL-ETTCN-65	6.5	10pc/box, 100pcs/CTN
YL-ETTCN-35	3.5	10pc/box, 100pcs/CTN	YL-ETTCN-70	7.0	10pc/box, 100pcs/CTN
YL-ETTCN-40	4.0	10pc/box, 100pcs/CTN	YL-ETTCN-75	7.5	10pc/box, 100pcs/CTN
YL-ETTCN-45	4.5	10pc/box, 100pcs/CTN	YL-ETTCN-80	8.0	10pc/box, 100pcs/CTN
YL-ETTCN-50	5.0	10pc/box, 100pcs/CTN	YL-ETTCN-85	8.5	10pc/box, 100pcs/CTN
YL-ETTCN-55	5.5	10pc/box, 100pcs/CTN	YL-ETTCN-90	9.0	10pc/box, 100pcs/CTN
YL-ETTCN-60	6.0	10pc/box, 100pcs/CTN	YL-ETTCN-95	9.5	10pc/box, 100pcs/CTN

REINFORCED ENDOTRACHEAL TUBE CUFFED



- Features :
- 1. There is a metal spring embedded inside the tube wall to avoid kinking.
  - 2. Balloon is made from famous COLORITE brand PVC material (DEHP Free available)
  - 3. X ray radio opaque line through
  - 4. Evenly balloon shape
  - 5. With Balloon

07

Product NO	Size	ID(mm)	suitable for Age	Product NO	Size	ID(mm)	suitable for Age
YL-ETTCR-30	#3.0	3.0	Pediatric weight 3kg	YL-ETTCR-65	#6.5	6.5	8-14 years
YL-ETTCR-35	#3.5	3.5	Pediatric 0-6 months	YL-ETTCR-70	#7.0	7.0	Female adult
YL-ETTCR-40	#4.0	4.0	6-12months	YL-ETTCR-75	#7.5	7.5	Female adult
YL-ETTCR-45	#4.5	4.5	1-3 years	YL-ETTCR-80	#8.0	8.0	Male adult
YL-ETTCR-50	#5.0	5.0	3-8 years	YL-ETTCR-85	#8.5	8.5	Male adult
YL-ETTCR-55	#5.5	5.5	3-8 years	YL-ETTCR-90	#9.0	9.0	Male adult
YL-ETTCR-60	#6.0	6.0	8-14 years	YL-ETTCR-100	#10.0	10.0	Male adult or Animal

REINFORCED ENDOTRACHEAL TUBE UNCUFFED

SKULL-ETUR



Brief Introduction:

- 1. Clear printing and opaque line through the tube
- 2. Smooth round tip minimize stimulation to patients during intubation
- 3. Standard 15mm semi-stick connector
- 4. Hard blister peelable packaging let E.O. Gas release easily
- 5. Complete specifications for different ages patients

Product NO	Size	ID(mm)	suitable for Age	Product NO	Size	ID(mm)	suitable for Age
YL-ETTUR-30	#3.0	3.0	Pediatric weight 3kg	YL-ETTUR-65	#6.5	6.5	8-14 years
YL-ETTUR-35	#3.5	3.5	Pediatric 0-6 months	YL-ETTUR-70	#7.0	7.0	Female adult
YL-ETTUR-40	#4.0	4.0	6-12months	YL-ETTUR-75	#7.5	7.5	Female adult
YL-ETTUR-45	#4.5	4.5	1-3 years	YL-ETTUR-80	#8.0	8.0	Male adult
YL-ETTUR-50	#5.0	5.0	3-8 years	YL-ETTUR-85	#8.5	8.5	Male adult
YL-ETTUR-55	#5.5	5.5	3-8 years	YL-ETTUR-90	#9.0	9.0	Male adult
YL-ETTUR-60	#6.0	6.0	8-14 years	YL-ETTUR-100	#10.0	10.0	Male adult or Animal

ENDOTRACHEAL TUBE WITH SUCTION LUMEN



Usage: A tracheal tube is a catheter that is inserted into the tracheal for the primary purpose of establishing and maintaining a patent airway and to ensure the adequate exchange of oxygen and carbon dioxide. ET tube, ETT tube are usually used as short name. Suction lumen is used for sputum suction.

Product features and advantages:

- 1. Material: balloon is made from famous COLORITE brand PVC material to ensure high quality
- 2. Symmetrical balloon shape
- 3. Clear printing and Opaque line through the tube
- 4. Smooth murphy stylets at the end
- 5. Smooth round tip minimize stimulation to patients during intubation
- 6. Standard 15mm semi-stick connector

Product NO	Cuff Size	Box/CTN	Product NO	Cuff Size	Box/CTN
YL-ETTCS-50	5.0	10pc/box, 100pcs/CTN	YL-ETTCS-80	8.0	10pc/box, 100pcs/CTN
YL-ETTCS-55	5.5	10pc/box, 100pcs/CTN	YL-ETTCS-85	8.5	10pc/box, 100pcs/CTN
YL-ETTCS-60	6.0	10pc/box, 100pcs/CTN	YL-ETTCS-90	9.0	10pc/box, 100pcs/CTN
YL-ETTCS-65	6.5	10pc/box, 100pcs/CTN	YL-ETTCS-95	9.5	10pc/box, 100pcs/CTN
YL-ETTCS-70	7.0	10pc/box, 100pcs/CTN	YL-ETTCS-100	10.0	10pc/box, 100pcs/CTN
YL-ETTCS-75	7.5	10pc/box, 100pcs/CTN			

SILICONE LARYNGEAL MASK



Features and benefits:

- 1. Made of pure medical grade silicone, have good biocompatibility, non-toxic
- 2. Exclusive soft seal cuff can be insert comfortable,minimizing potential trauma and increasing sealing
- 3. Strengthened neck and tip eases insertion and prevents folds during insertion
- 4. Have different sizes, suitable for infant, children and adult

Product NO	Size(Fr)	Maximum Cuff Inflation Volumes ( ml )	Patient Weight(KG)
YL-SLMA-10	1.0	4	0-5
YL-SLMA-15	1.5	7	5-10
YL-SLMA-20	2.0	10	10-20
YL-SLMA-25	2.5	14	20-30
YL-SLMA-30	3.0	20	30-50
YL-SLMA-40	4.0	30	50-70
YL-SLMA-50	5.0	40	70-100



REUSABLE / REINFORCED LARYNGEAL MASK



Features and benefits:

- 1. Made of pure medical grade silicone, have good biocompatibility, non-toxic
- 2. Exclusive soft seal cuff can be insert comfortable, minimizing potential trauma and increasing sealing
- 3. Strengthen neck and tip eases insertion and prevents folds during insertion
- 4. Reinforced LMA specially designed for ENT, opthal mic, dental and other head and neck surgeries
- 5. Kink-free wire reinforced tube eliminates the risk of airway tube occlusion
- 3. Reusable types must be cleaning and disinfection or sterilization before use, use only 40 times
- 7. Every product with unique serial number and record card to avoid confusion and ensures the correct record offer use
- 8. Have different sizes, suitable for infant, children and adult.

Product NO	Size(Fr)	Maximum Cuff Inflation Volumes ( ml )	Patient Weight(KG)
YL-RLMA-10	1.0	4	0-5
YL-RLMA-15	1.5	7	5-10
YL-RLMA-20	2.0	10	10-20
YL-RLMA-25	2.5	14	20-30
YL-RLMA-30	3.0	20	30-50
YL-RLMA-40	4.0	30	50-70
YL-RLMA-50	5.0	40	70-100

PVC LARYNGEAL MASK



Specifications:

- 1. The tube and cuff are made of medical-grade PVC
- 2. Specially designed shape coincides with the laryngopharynx well, reducing stimulation to patient body and improving the cuff seal.
- 3. Convenient and safe
- 4. Sterilized, single use only
- 5. Seven sizes, suitable for adult, children

12

Product NO	Size(Fr)	Maximum Cuff Inflation Volumes ( ml )	Patient Weight(KG)
YL-PVC LMA-10	1.0	4	0-5
YL-PVC LMA-15	1.5	7	5-10
YL-PVC LMA-20	2.0	10	10-20
YL-PVC LMA-25	2.5	14	20-30
YL-PVC LMA-30	3.0	20	30-50
YL-PVC LMA-40	4.0	30	50-70
YL-PVC LMA-50	5.0	40	70-100

ENHANCED SUCTION TUBE



- Products Description
- Medical Closed Suction System Catheter 72 Hours
1. The closed suction tube has designed for proven effective in preventing infections, reducing crosscontamination, reducing intensive care unit days and patient costs.
  2. Providing Quality Solutions for respiratory care.
  3. Sterile,individual PU protective sleeve of closed suction system can protect the caregivers from cross infection.
  4. With isolation valve for effective VAP control.
  5. Individually wrapped to stay fresh.
  6. Closed suction tube sterilized by EO gas,latex free and for single-use.
  7. Double swivel connectors reduce strain on ventiltor tubing.

Product NO	Size	Type	Length for option	Attachment for option
YL-CSC-06F-72H	Fr06	72 hour	36cm, 54cm, 60cm	Y connector
YL-CSC-08F-72H	Fr08		36cm, 54cm, 60cm	Y connector
YL-CSC-10F-72H	Fr10		36cm, 54cm, 60cm	Y connector
YL-CSC-12F-72H	Fr12		54cm, 60cm	/
YL-CSC-14F-72H	Fr14		54cm, 60cm	/
YL-CSC-16F-72H	Fr16		54cm, 60cm	/
YL-CSC-18F-72H	Fr18		54cm, 60cm	/

DISPOSABLE CLOSED SUCTION CATHETER

14



- Products Description
- Medical Closed Suction System Catheter 24 Hours
1. Raditional secretion suction should cut off patient’s breathing airway, so nurse has to finish suction in 30 seconds.
  2. Suppose closed suction system is reused for 48 hours, the average cost is almost the same as using traditional suction catheters.
  3. It reduce cross-infection during extraction secretion

Product NO	Size	Type	Length for option	Attachment for option
YL-CSC-06F-24H	Fr06	24 hour	36cm, 54cm, 60cm	Y connector
YL-CSC-08F-24H	Fr08		36cm, 54cm, 60cm	Y connector
YL-CSC-10F-24H	Fr10		36cm, 54cm, 60cm	Y connector
YL-CSC-12F-24H	Fr12		54cm, 60cm	/
YL-CSC-14F-24H	Fr14		54cm, 60cm	/
YL-CSC-16F-24H	Fr16		54cm, 60cm	/
YL-CSC-18F-24H	Fr18		54cm, 60cm	/

PVC SUCTION CATHETER



Usage: The suction catheter and suction tubing are connected to a collection container (suction canister) and a device that generates suction. This is a single use catheter which is used most often when the patient does not have an artificial airway such as an endotracheal tube or tracheost

Features:

- 1. Round and polished catheter tip
- 2. Tube hardness is changeable according to patient age.
- 3. Color-coded connector at proximal end for instant identification of catheter size
- 4. Available with thumb control device for positive control during suction

15

Product NO	Size	O.D.	Color	Length
YL-SC05G	Fr5	1.7mm	Gray	380 ~ 600 mm
YL-SC06G	Fr6	2.0mm	Light green	380 ~ 600 mm
YL-SC08G	Fr8	2.7mm	Light blue	400 ~ 600 mm
YL-SC10G	Fr10	3.3mm	Black	400 ~ 600 mm
YL-SC12G	Fr12	4.0mm	White	400 ~ 600 mm
YL-SC14G	Fr14	4.7mm	Green	400 ~ 600 mm
YL-SC16G	Fr16	5.3mm	Orange	400 ~ 600 mm
YL-SC18G	Fr18	6.0mm	Red	400 ~ 600 mm
YL-SC20G	Fr20	6.7mm	Yellow	400 ~ 600 mm
YL-SC22G	Fr22	7.3mm	Purple	400 ~ 600 mm



SUCTION CONNECTING TUBE + YANKAUER ( SUCTION SET )

Features:

- 1. Made from transparent material for better visualization
- 2. Striated wall of tube provide superior strength and anti-kinking
- 3. Multiple choices of length
- 4. Available with small connector which can connect with suction catheter
- 5. Flared connector could connect with smooth suction yankauer handle
- 6. Yankauer suction handle with or without vacuum control



16

Product NO	Depth	Product NO	Depth
YL-SSYT180	1.8m	YL-SSYT300	3.0m
YL-SSYT200	2.0m	YL-SSYT360	3.6m
YL-SSYT210	2.1m	YL-SSYT500	5.0m
YL-SSYT250	2.5m		

PVC RYLE TUBE / NASOGASTRIC TUBE



Features:

- 1. With distal end open tip ( Closed tip available), atraumatic
- 2. Ticker tube can be also used to suck out stomach fluid for testing
- 3. X – ray line detectable
- 4. Kink resistant
- 5. Sterile

Product NO	size	color	Length
YL-RTF10	Fr10	Black	120-140cm
YL-RTF12	Fr12	White	120-140cm
YL-RTF14	Fr14	Green	120-140cm
YL-RTF16	Fr16	Orange	120-140cm
YL-RTF18	Fr18	Red	120-140cm
YL-RTF20	Fr20	Yellow	120-140cm

TRACHEOSTOMY TUBE WITH CUFF



- Brief Introduction:
- 1. Material: Balloon is made from famous COLORITE brand PVC material to ensure high quality
  - 2. Symmetrical balloon shape
  - 3. Clear printing and opaque line through the tube
  - 4. Smooth round tip minimize stimulation to patients during intubation
  - 5. Standard 15mm semi-stick connector
  - 6. Hard blister peelable packaging let E.O. Gas release easily
  - 7. Complete specifications for different ages patients

Product NO	Cuff Size	Box/CTN	Product Number	Cuff Size	Box/CTN
YL-TOTC-50	5.0	10pc/box, 100pcs/CTN	YL-TOTC-70	7.0	10pc/box, 100pcs/CTN
YL-TOTC-55	5.5	10pc/box, 100pcs/CTN	YL-TOTC-75	7.5	10pc/box, 100pcs/CTN
YL-TOTC-60	6.0	10pc/box, 100pcs/CTN	YL-TOTC-80	8.0	10pc/box, 100pcs/CTN
YL-TOTC-65	6.5	10pc/box, 100pcs/CTN	YL-TOTC-90	9.0	10pc/box, 100pcs/CTN

TRACHEOSTOMY TUBE WITHOUT CUFF

Brief Introduction:

- 1. Clear printing and opaque line through the tube
- 2. Smooth round tip minimize stimulation to patients during intubation
- 3. Standard 15mm semi-stick connector
- 4. Hard blister peelable packaging let E.O. Gas release easily
- 5. Complete specifications for different ages patients

19

Product NO	Cuff Size	Box/CTN	Product Number	Cuff Size	Box/CTN
YL-TOTU-30	3.0	10pc/box, 100pcs/CTN	YL-TOTU-65	6.5	10pc/box, 100pcs/CTN
YL-TOTU-35	3.5	10pc/box, 100pcs/CTN	YL-TOTU-70	7.0	10pc/box, 100pcs/CTN
YL-TOTU-40	4.0	10pc/box, 100pcs/CTN	YL-TOTU-75	7.5	10pc/box, 100pcs/CTN
YL-TOTU-45	4.5	10pc/box, 100pcs/CTN	YL-TOTU-80	8.0	10pc/box, 100pcs/CTN
YL-TOTU-50	5.0	10pc/box, 100pcs/CTN	YL-TOTU-85	8.5	10pc/box, 100pcs/CTN
YL-TOTU-55	5.5	10pc/box, 100pcs/CTN	YL-TOTU-90	9.0	10pc/box, 100pcs/CTN
YL-TOTU-60	6.0	10pc/box, 100pcs/CTN			

NEBULIZER OXYGEN MASK

Features:

- 1、 Made of medical grade PVC
- 2、 With adjustable nose clip
- 3、 With elastic strap
- 4、 With 2m tubing
- 5、 With 6ml Nebulizer cup
- 6、 Size: S, M, L, XL

Product NO	Size	Cs(pcs)
YL-NOM-XL	Adult Elongated/XL	100
YL-NOM-L	Adult Standard/L	100
YL-NOM-M	Pediatric Elongated/M	100
YL-NOM-S	Pediatric Standard/S	100



DISPOSABLE OXYGEN MASK

Features:

- 1. Brand Name: EasyThru or OEM, ODM
- 2. Certification: CE ISO13485
- 3. Model Number: S, M, L, XL
- 4. Material: Medical Grade PVC
- 5. Size: S, M, L, XL Suitable For Pediatric And Adult
- 6. Tube Length: 2m, 2.1m

Product NO	Size	Cs(pcs)
YL-OM-XL	Adult Elongated/XL	100
YL-OM-L	Adult Standard/L	100
YL-OM-M	Pediatric Elongated/M	100
YL-OM-S	Pediatric Standard/S	100





TRACHEOSTOMY MASK

Features Of Tracheostomy Mask

- 1. Be used to deliver oxygen gas to tracheostomy patients;
- 2. Be worn around patient's neck over the tracheostomy tube.
- 3. PE packing with label inside;
- 4. Tube connector swivels 360 degree for different position of patients.
- 5. Both adult size and pediatric size available;
- 6. CE, ISO approved.

(ADULT TYPE)



(CHILD TYPE)



Product No	Size(Fr)	Msak	Cs(pcs)
YL-TMA	Adult	Without Tubing	100
YL-TMC	Child	Without Tubing	100

2 WAY LATEX FOLEY CATHETER WITH SILICONE COATED



- Features:
- 1. All products pass cytotoxicity testing according to ISO 10993
  - 2. Latex material 100% imported from Malaysia, specially used to produce latex Foley catheter
  - 3. Unique designed rubber valve ensure leakage rate 40%-60% less than others in the market

22

Product NO	Size(Fr)	Valve	Balloon Volume( ml )	Length(mm)
YL-LFC2-06	6	Rubber/plastic	3ml / 3-5ml	270-280
YL-LFC2-08	8	Rubber/plastic	3-5ml	270-280
YL-LFC2-10	10	Rubber/plastic	3-5ml / 5ml	270-280
YL-LFC2-12	12	Rubber/plastic	5ml/ 10ml/ 5-10ml/ 5-15ml	390-400
YL-LFC2-14	14	Rubber/plastic	5ml/ 10ml/ 5-10ml/ 5-15ml	
YL-LFC2-16	16	Rubber/plastic	5-10ml / 5-15ml / 30ml / 30-50ml	
YL-LFC2-18	18	Rubber/plastic	5-10ml / 5-15ml / 30ml / 30-50ml	
YL-LFC2-20	20	Rubber/plastic	5-15ml / 30ml / 30-50ml	
YL-LFC2-22	22	Rubber/plastic	5-15ml / 30ml / 30-50ml	
YL-LFC2-24	24	Rubber/plastic	30ml / 30-50ml	
YL-LFC2-26	26	Rubber/plastic	30ml / 30-50ml	

3 WAY LATEX FOLEY CATHETER WITH SILICONE COATED



- Features:
- The third way is used for irrigation or medicine injection
- 1. All products pass cytotoxicity testing according to ISO 10993
  - 2. Latex material 100% imported from Malaysia, specially used to produce latex Foley catheter
  - 3. Unique designed rubber valve ensure leakage rate 40%-60% less than others in the market

Product NO	Size (Fr)	Valve	Balloon Volume(ML)	Length (mm)
YL-LFC3-16	16	Rubber/plastic	30ml/30-50ml	390-400
YL-LFC3-18	18	Rubber/plastic	30ml/30-50ml	
YL-LFC3-20	20	Rubber/plastic	30ml/30-50ml	
YL-LFC3-22	22	Rubber/plastic	30ml/30-50ml	
YL-LFC3-24	24	Rubber/plastic	30ml/30-50ml	
YL-LFC3-26	26	Rubber/plastic	30ml/30-50ml	

## HAEMOSTATIC CATHETER (COUVELAIRE TIP DUFOR TIP , MERCIER TIP)

Usage: Prostatectomy, Prostate cancer surgery, Bladder cancer resection and urinary surgery operation

Advantages:

1. After operation, this catheter is user for Quick compr-ession hemostasis by big balloon, Quick suction away blood clot, Drainage hard substance.Irrigation bladder, urine drainage, Drug Injection.
2. Compression hemos tasis: balloon injextion volume up to 80ml Balloon blasting limit to 250ml




## ENHANCED HAEMOSTATIC CATHETER (COUVELAIRE TIP DUFOR TIP, MERCIER TIP)

Usage: Prostatectomy, Prostate cancer surgery, Bladder cancer resection and urinary surgery operation

Advantages:

1. After operation, this catheter is user for Quick compression hemostasis by big balloon, Quick suction away blood clot, Drainage hard substance.Irrigation bladder, urine drainage, Drug injection.
2. Compression hemos tasis: balloon injextion volume up to 80ml Balloon blasting limit to 250ml
3. Medical grade spring through tube wall, anti-collapse during clot suction under negative pressure 80mpa.



DOUBLE BALLOON 3 WAY LATEX URINARY CATHETERS



Usage:  
Front balloon fasten placement, the back one 50ml is used to stop blood by pressure

Product NO	Size (FR/CH)	Color Code	O.D. (mm)	Balloon Size
YL-DB3-18	Fr18	Red	0.6mm	30ml + 50ml
YL-DB3-20	Fr20	Yellow	6.6mm	30ml + 50ml
YL-DB3-22	Fr22	Purple	7.3mm	30ml + 50ml
YL-DB3-24	Fr24	Blue	8.0mm	30ml + 50ml
YL-DB3-26	Fr26	Pink	8.7mm	30ml + 50ml

2 WAY TIEMMANN TIP LATEX FOLEY CATHETER



Features:

- 1. All products pass cytotoxicity testing according to ISO 10993.
- 2. Latex material 100% imported from Malaysia, specially used to produce latex foley catheter
- 3. Unique design rubber valve ensure leakage rate 40%-60% less than others in the market
- 4. Stable flow in both drainage lumen and sterilized water injection lumen
- 5. Symmetry balloon after inflating and No burst balloon according to testing
- 6. Smooth and symmetry murphy eyelets
- 7. Certificates and tests: ISO, CE, Free Sales Certificate, Product Registration, Cytotoxicity testing
- 8. Hydrophilic coating is available

Product NO	Size	Color Code	O.D. (mm)	Balloon Size
YL-TTLFC2-12	Fr12	White	4.0	5ml - 10ml, 30ml, 30ml - 50ml
YL-TTLFC2-14	Fr14	Green	4.6	5ml - 10ml, 30ml, 30ml - 50ml
YL-TTLFC2-16	Fr16	Orange	5.3	5ml - 10ml, 30ml, 30ml - 50ml
YL-TTLFC2-18	Fr18	Red	6.0	5ml - 10ml, 30ml, 30ml - 50ml
YL-TTLFC2-20	Fr20	Yellow	6.7	5ml - 10ml, 30ml, 30ml - 50ml
YL-TTLFC2-22	Fr22	Purple	7.3	5ml - 10ml, 30ml, 30ml - 50ml
YL-TTLFC2-24	Fr24	Blue	8.0	5ml - 10ml, 30ml, 30ml - 50ml

2 WAY SILICONE FOLEY CATHETER



Features:

- 1. Made of 100% imported medical-grade silicone
- 2. The retention time is within 29 days
- 3. Radio opaque line through the length for x-ray visualization
- 4. Soft and uniformly inflated balloon makes tube sit well against the bladder.
- 5. Color-coded check valve for identification of different sizes.
- 6. Length of the Foley catheter: Children: 310mm (with guide wire), Adult: 410mm.

Product NO	Size(Fr)	Balloon volume(ML)	Length
YL-SFC2-06	6	3ml	310mm
YL-SFC2-08	8	3-5ml	
YL-SFC2-10	10	3-5ml	
YL-SFC2-12	12	5ml/5-10ml/10ml/30ml	410mm
YL-SFC2-14	14		
YL-SFC2-16	16		
YL-SFC2-18	18		
YL-SFC2-20	20		
YL-SFC2-22	22		
YL-SFC2-24	24		

3 WAY SILICONE FOLEY CATHETER

28



- Features:
- 1. Made of 100% imported medical-grade silicone
  - 2. The retention time is within 29 days
  - 3. Radio opaque line through the length for x-ray visualization
  - 4. Soft and uniformly inflated balloon makes tube sit well against the bladder.
  - 5. Color-coded check valve for identification of different sizes.
  - 6. Length of the Foley catheter: 400-410mm

Product NO	Size(Fr)	Balloon volume(ml)	Length (mm)
YL-SFC3-14	14	30ml	410 mm
YL-SFC3-16	16		
YL-SFC3-18	18		
YL-SFC3-20	20		
YL-SFC3-22	22		
YL-SFC3-24	24		



PVC INTERMITTENT URETHRAL CATHETER  
PVC NELATON CATHETER




- Features:
- 1. Medical grade PVC, non- toxic
  - 2. Smooth, closed tip with two lateral eyelets
  - 3. Color coded connector easy to identify size
  - 4. Used for short time intermittent urine drain
  - 5. Sterile, single use
  - 6. DEHP free for option

Product NO	Type	Size	Length	Individual packaging	Carton
YL-PVC NC-M	Male	Fr6-Fr22	40cm	PE	100pcs/box 600pcs/CTN
				Blister	
YL-PVC NC-F	Female	Fr6-Fr22	20cm	PE	100pcs/box 600pcs/CTN
				Blister	
YL-PVC NC-T	Tiemann	Fr6-Fr22	40cm	PE	100pcs/box 600pcs/CTN
				Blister	

FIRE-POLISHED INTERMITTENT URETHRAL CATHETER

30



Features:

- 1. Fire – polished eyelets are much smoother and comfortable than regular eyelets
- 2. DEHP Free, Latex Free
- 3. Smooth, closed tip
- 4. Color coded connector easy to identify size
- 5. Used for short time intermittent urine drain
- 6. Sterile, single use

Product NO	Type	Size	Length	Individual packaging	Carton
YL-FIUC-M	Male	Fr6-Fr22	40cm	PE	100pcs/box 600pcs/CTN
				Blister	
YL-FIUC-F	Female	Fr6-Fr22	20cm	PE	100pcs/box 600pcs/CTN
				Blister	
YL-FIUC-T	Thiemann	Fr6-Fr22	40cm	PE	100pcs/box 600pcs/CTN
				Blister	

PLASTIBELL CIRCUMCISION



Usage: Plastibell circumcision is mainly used in urology surgery, intended for circumcision of patients with redundant prepuce and phimosis. It makes the complex circumcision surgery become simple and convenient, fast and safe.

Size: 11/12/13/14/15/16/17/18/18A/19

Advantages:

- 1. Safety and Accurate, single use only avoid cross infection.
- 2. Cut, blend automatically
- 3. Simple operate, finished in 3-5 minutes
- 4. No Suturing, Minimal Bleeding, Little Pain
- 5. Rapid Recovery
- 6. ISO, CE certificate

Product NO	Specification	Size	Comosition:
YL-PC-11	11#	11	Stemloop(withhandle), suture string
YL-PC-12	12#	12	
YL-PC-13	13#	13	
YL-PC-14	14#	14	
YL-PC-15	15#	15	
YL-PC-16	16#	16	
YL-PC-17	17#	17	
YL-PC-18	18#	18	
YL-PC-18A	18A#	18A	
YL-PC-19	19#	19	

CIRCUMCISION DEVICE



**Usage:** In circumcision, the foreskin is removed to exposed the head of the penis. It's a quick procedure that causes very little bleeding and stitches aren't needed.

- Advantages:**
- 1. Safety and Accurate, single use only avoid cross infection.
  - 2. Cut, blend automatically
  - 3. Simple operate, finished in 5-10 minutes
  - 5. Full sizes for all ages
  - 6. No Suturing, Minimal Bleeding, Little Pain
  - 7. Rapid Recovery
  - 8. ISO, CE certificate

Product NO	Diameter	Number of staples ( pieces)	Staples Material
YL-CD-12	12mm	9	Stainless Steel or Titanium alloy
YL-CD-15	15mm	12	
YL-CD-18	18mm	12	
YL-CD-21	21mm	15	
YL-CD-26	26mm	15	
YL-CD-30	30mm	20	
YL-CD-32	32mm	22	
YL-CD-36	36mm	22	

URINE DRAINAGE BAGS



Simple urine bag  
Pull – push valve



Simple urine bag  
Cross valve



Luxury urine bag  
Air-Vent and  
Anti-Reflux Valve



Precision urine bag  
T-Tap water outlet valve  
urine meter with twist  
valveanti-reflux valve

Product NO	Name	Capacity	Valve type
YL-SUB-P	Simple urine bag	1000ml, 2000ml	Pull – push valve
YL-SUB-C	Simple urine bag	1000ml, 2000ml	Cross valve
YL-LUB	Luxury urine bag	1500 ml, 2000ml	Air-Vent and Anti-Reflux Valve
YL-PUB	Precision urine bag	2000ml bag+200ml bottle	T-Tap water outlet valve urine meter with twist valve anti-reflux valve
		2500ml bag+500ml bottle	

MALE EXTERNAL CATHETER

34



Latex External Catheter

Product NO	Size	Type
YL-LEC-S	20mm	S
YL-LEC-M	25mm	M
YL-LEC-L	30mm	L
YL-LEC-XL	35mm	XL

Silicone External Catheter

Product NO	Size	Type
YL-SEC-S	25mm	S
YL-SEC-M	30mm	M
YL-SEC-L	35mm	L
YL-SEC-XL	40mm	XL



INSULIN PEN NEEDLE



Usage: Insulin pen needle is together used with insulin pen, which offer several significant advantages over insulin syringes, such as easy of handling, accuracy and comfortable.

Features and Advantages:

- 1. Comfortable and painless.
- 2. Three bevel needle for a gentle and minimum puncture.
- 3. Silicone lubricant coating on insulin pen needles for maximum injection comfort.
- 4. Safe and convenient to carry
- 5. Non-Toxic, non- pyrogen needle guaranteed safety
- 6. Easy to use by simply twisting on & off
- 7. Full size for all people
- 8. Fit major brands injection pens

Product NO	Needle Size	Needle Length	Count/Box
YL-IPN29G12	29G 0.33mm	12mm	100pc/box,5000pcs/CTN
YL-IPN30G10	30G 0.30mm	10mm	100pc/box,5000pcs/CTN
YL-IPN30G08	30G 0.30mm	8mm	100pc/box,5000pcs/CTN
YL-IPN31G08	31G 0.25mm	8mm	100pc/box,5000pcs/CTN
YL-IPN31G06	31G 0.25mm	6mm	100pc/box,5000pcs/CTN
YL-IPN31G05	31G 0.25mm	5mm	100pc/box,5000pcs/CTN
YL-IPN32G04	32G 0.23mm	4mm	100pc/box,5000pcs/CTN
YL-IPN33G04	32G 0.20mm	4mm	100pc/box,5000pcs/CTN

INSULIN SYRINGE



- 1. Insulin syringes are made specifically for self injections and have friendly features: Shorter needles, as infulin injections are subcutaneous (under the skin) rather than intramuscular
- 2. Finer gauge needles for less pain
- 3. Markings in insulin units to simplify drawing a measured dose of insulin
- 4. Easy to read the scale on syringe
- 5. Clear Barrel make sure the air bubble are easy to be idendtified
- 6. mooth action of the plunger when drawing up insulin
- 7. Sterilized, single use only, safety

Product NO	Insulin Syringe	Syringe Size	Needle Size	Product NO	Insulin Syringe	Syringe Size	Needle Size
YL-IS100-032912	U - 100	0.3ml	29*1/2''	YL-IS40-052912	U - 40	0.5ml	29*1/2''
YL-IS100-052912		0.5ml	29*1/2''	YL-IS40-053012		0.5ml	30*1/2''
YL-IS100-053012		0.5ml	30*1/2''	YL-IS40-0530516		0.5ml	30*5/16''
YL-IS100-0530516		0.5ml	30*5/16''	YL-IS40-102912		1.0ml	29*1/2''
YL-IS100-0531516		0.5ml	31*5/16''	YL-IS40-103012		1.0ml	30*1/2''
YL-IS100-102912		1.0ml	29*1/2''	YL-IS40-1030516		1.0ml	30*5/16''
YL-IS100-103012		1.0ml	30*1/2''	YL-IS40-103112		1.0ml	31*1/2''
YL-IS100-1030516		1.0ml	30*5/16''	YL-IS40-1031516		1.0ml	31*5/16''
YL-IS100-103112		1.0ml	31*1/2''				
YL-IS100-1031516		1.0ml	31*5/16''				



SAFETY BLOOD LANCET



Features and advantages:

- 1. User safety: needle is fully shielded before and after use.
- 2. Patient safety: Sterilized, needle protected by safety cap, built-in needle protection system.
- 3. Easy handling: Only 2-3 steps finish blood sampling
- 4. Patient comfort: Sharp needle, little pain feeling
- 5. Application: blood glucose and other in diagnostic tests.
- 6. Protection: protect the healthcare professionals from the risk of HBV, HCV, HIV.

Product NO	Color	Size	Depth
YL-SBL-21	Yellow	21G	2.2mm
YL-SBL-23	Green	23G	2.2mm
YL-SBL-26	Blue	26G	1.8mm
YL-SBL-28	Purple	28G	1.8mm
YL-SBL-30	Pink	30G	1.5mm

NEEDLE HOLDER



Safety cover lancet holder

Advantage of our needle holders

- 1. Transparent plastic material is easy for doctors to observe blood collection status of vacuum tubes.
- 2. Double screwing juncture can fasten blood collection needle tightly
- 3. Green cover offer maximal protection
- 4. General luer lock, fit for all kinds pen like needles

Snap-in lancet frame

Advantage of our needle holders

- 1. Transparent plastic material is easy for doctors to observe blood collection status of vacuum tubes.
- 2. Double screwing juncture can fasten blood collection needle tightly.
- 3. Green cover offer maximal protection
- 4. General luer lock, fit for all kinds pen like needles



Product NO	Name	Material	Model Number
YL-CBNH-SCLH	Safety cover lancet holder	PP	General Type
YL-CBNH-SLF	Snap-in lancet frame	PP	General Type / Snap-type

SILICONE STOMACH TUBE



Stomach tube/ Feeding tube is used to deliver food, nutrients, medication, or other material into the stomach, or drain undesirable contents from the stomach, or decompress the stomach. It is inserted through patient's mouth or nose.

- Features:
- 1. Made from silicone
  - 2. X-ray through the tube length
  - 3. Clear depth printing
  - 4. E.O. sterilized, single use only

Product NO	Size (Fr/Ch)	Length(cm)	Marked Position(cm)
YL-SST-06	6	120	40,45,50,55,60,70,80
YL-SST-08	8		
YL-SST-10	10		
YL-SST-12	12		
YL-SST-14	14		
YL-SST-16	16		
YL-SST-18	18		
YL-SST-20	20		

TPU FEEDING TUBE

40




Features:

- 1. Manufactured from 100% medical grade TPU for superior biocompatibility
- 2. Proper softness and excellent flexibility for enhanced patient comfort
- 3. Tube with radio-opaque line throughout the length
- 4. Distal end is coned with 4 lateral eyes provide efficient drainage and irrigation and minimize blockage of the lumen
- 5. Depth markings in 10cm intervals for accurate placement
- 6. Universal funnel shaped connector with attached cap for secure sealing
- 7. Supplied in peelable pouch for aseptic handling

Product NO	Size (Fr/Ch)	Length(cm)	Marked Position(cm)
YL-PUFT-06	6	120	40,50,60,70,80
YL-PUFT-08	8		
YL-PUFT-10	10		
YL-PUFT-12	12		
YL-PUFT-14	14		
YL-PUFT-16	16		
YL-PUFT-18	18		
YL-PUFT-20	20		

PVC FEEDING TUBE



Features:

- 1、 Made from medical grade PVC
- 2、 DEHP FREE available
- 3、 line available
- 4、 kink resistance
- 5、 Smooth round tip
- 6、 Cover connector
- 7、 Graduated scale printing
- 8、 Frosted or Transparent tube for option
- 9、 Color coded connector

Product NO	Size (Fr/Ch)	Length(cm)	Product NO	Size (Fr/Ch)	Length(cm)
YL-PFT-05	5	50	YL-PFT-14	14	50
YL-PFT-06	6		YL-PFT-16	16	
YL-PFT-08	8		YL-PFT-18	18	
YL-PFT-10	10		YL-PFT-20	20	
YL-PFT-12	12				

N95 DUSTY MASK

42



Type: Unfoldable  
Vlave: with valve or without valve  
Elastic ear loop  
Adjustable nose clip




Type: Foldable  
Vlave: with valve or without valve  
Elastic ear loop  
Adjustable nose clip

Product NO	Type	valve
YL-N95DM-U	Unfoldable	No valve
YL-N95DM-F	Foldable	With valve

DISPOSABLE MEDICAL FACE MASK (TYPE IIR, EARLOOP) SURGICAL FACE MASK

Figure 1



Usage: Type IIR disposable Medical Face Mask also named surgical mask, which are intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements, to protect the wearer from infective agents and, additionally, in certain circumstances to protect the wearer against splashes of potentially contaminated liquids. Surgical Face Mask also prevent the spread of dander and respiratory microorganisms to patient’s open surgical wound. Thus it playing a role of two-way biological protection.

Features:

- 1. 3 plys, size 17.5\*9.5cm
- 2. BFE ≥98%, Differential pressure< 40, Splash resistance pressure (kpa)≥16.0; Microbial cleanliness ( cfu/g)≤30 ;
- 3. Conform to EN14683:2019

Product NO	Specification		Count/Box	Count/ CTN	CTN SIZE
YL-DFM-IIR-LOP01	<ul style="list-style-type: none"><li>• 3 plys</li><li>• Size 17.5*9.5cm</li><li>• BFE≥98%</li><li>• Differential pressure&lt; 40</li></ul>	Non Sterile	50 pcs / box	2000 pcs/CTN	51*38*42cm
YL-DFM-IIR-LOP02	<ul style="list-style-type: none"><li>• Splash resistance pressure ( kpa)≥16.0;</li><li>• Microbial cleanliness ( cfu/g)≤30;</li><li>• EN14683:2019</li></ul>	Sterile	50 pcs/ box	2000 pcs/CTN	51*38*42cm

DISPOSABLE MEDICAL FACE MASK (TYPE IIR, TIE) SURGICAL FACE MASK



Usage: Type IIR disposable Medical Face Mask also named surgical mask, which are intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements, to protect the wearer from infective agents and, additionally, in certain circumstances to protect the wearer against splashes of potentially contaminated liquids. Surgical Face Mask also prevent the spread of dander and respiratory microorganisms to patient's open surgical wound. Thus it playing a role of two-way biological protection.

Features:

- 1. 3 plys, size 17.5\*9.5cm
- 2. BFE ≥98%, Differential pressure< 40, Splash resistance pressure (kpa)≥16.0; Microbial cleanliness (cfu/g)≤30 ;
- 3. Conform to EN14683:2019

Product NO	Specification		Count/Box	Count/ CTN	CTN SIZE
YL-DFM-IIR-TIE01	<ul style="list-style-type: none"><li>• 3 plys</li><li>• Size 17.5*9.5cm</li><li>• BFE≥98%</li><li>• Differential pressure&lt; 40</li></ul>	Non Sterile	50 pcs / box	2000 pcs/CTN	51*38*42cm
YL-DFM-IIR-TIE02	<ul style="list-style-type: none"><li>• Splash resistance pressure ( kpa)≥16.0;</li><li>• Microbial cleanliness ( cfu/g)≤30;</li><li>• EN14683:2019</li></ul>	Sterile	50 pcs/ box	2000 pcs/CTN	51*38*42cm



## DISPOSABLE SURGICAL GLOVES

Sterilized Medical Surgical Latex Gloves Powder or Powder free.

Size: #6.0, #6.5, #7.0, #7.5, #8.0, #8.5, #9.0



## DISPOSABLE LATEX SURGICAL GLOVES

Advantages:

1. Extra strength provides additional protection from surgical debris.
2. Fully anatomical design to reduce hand fatigue.
3. Softness provides superior comfort and natural fit.
4. Micro-roughened surface provides excellent wet and dry grip.
5. Easy donning and helps prevent rolling back.
6. High strength and elasticity

45

## DISPOSABLE LATEX EXAMINATION GLOVES

Advantages:

1. Extra strength provides additional protection from surgical debris.
2. Fully anatomical design to reduce hand fatigue.
3. Softness provides superior comfort and natural fit.
4. Micro-roughened surface provides excellent wet and dry grip.
5. Easy donning and helps prevent rolling back.
6. High strength and elasticity





## **NANCHANG YILI MEDICAL INSTRUMENT CO., LTD**

Office: 101# Building, Greenland International Exp Center, Jiulong Road 1177, Honggutan District, Nanchang, Jiangxi, China

Factory: Jinxian Medical Industrial Park, Nanchang, Jiangxi, China

Tel: +86-791-88252326 Fax: +86-791-88251326 Post Code: 330046

E-mail: [sales@yilimed.com](mailto:sales@yilimed.com)

[www.yili-medical.com](http://www.yili-medical.com)

